

UNO NARROW

The Uno Narrow replaceable ball anchor screw, which is purchased separately, is available in three varying heights and its ball is of standard dimensions, making it compatible with our entire line of ball caps.

length

Type

Ø3mm
screw type implants

10mm

MO2-10300



11.5mm

MO2-11300




13mm

MO2-13300



16mm

MO2-16300




Ball attachment anchor:



MB-NRO05
Height 0.5mm




MB-NRO15
Height 1.5mm



MB-NRO25
Height 2.5mm

UNO NARROW

Tools




 Not include in the advance kit (MK-0035, MK-EI35 MK-IR35) and can be purchase seperately*

Procedure

(Ø3mm Implant)

Drill Speed (RPM)	1200-1500	900-1200	700-900	15-20
Diameter	Ø 1.90	Ø 2	Ø 2.40	Ø 3



Optional for dense bone